



CANAVERAL MODULAR SCREWS



The **FloSpine Canaverall™ System** offers a unique approach for initial screw placement followed by head body attachment . This provides the following benefits:

- **Improved identification and targeting of structural landmarks**
- **Enhanced visibility during innerbody work**
- **Support for easier distraction**
- **Flexibility in head body selection**

We've Got
Your Back™

Modular

CANAVERAL



MODULAR SCREWS

PRECISION IMPLANTS


- ✓ Multiple Headbody Styles
- ✓ Cannulated Screws Available
- ✓ Designed to Integrate with Robotic Assisted placement
- ✓ Comprehensive Screw Offering, 4.5mm-10.5mm




STREAMLINED INSTRUMENTS



Place Screws



Distract and Perform Innerbody Work



Place HeadBody and Rod

Manufactured in United States of America
* Other sizes not shown available upon request
Caution: Federal law (USA) restricts this device to sale by or on the order of a physician
561-705-3080
www.flospine.com
XR-1309-PL-US FloSpine, LCC 11/2014